

# The Photoelectric Effect

Enrique Rivera Jr.  
Physics Undergraduate,  
The University of Texas at Austin

April 7, 2024

## Abstract

The photoelectric effect is the emission of electrons from a metal surface when irradiated with photons whose energies are greater than the work function of the metal. This experiment aims to measure Planck's constant ( $h$ ) and the work function ( $\phi$ ) of the metal by analyzing the photoelectric effect. The experimental setup includes a mercury lamp as the photon source, a lens to focus the light, and an RCA 935 vacuum phototube. Interference filters were used to select photon wavelengths, and neutral density filters tested the dependence of photoelectron energy on light intensity. The results of the experiment were consistent with the theoretical predictions, and the values of Planck's constant and the work function were determined to be  $h = 8.55 \times 10^{-31} J \cdot s$  and  $\phi = 1.97 \times 10^{-22} J$ , respectively. Which come close to the accepted values of  $h = 6.63 \times 10^{-34} J \cdot s$  and  $\phi = 3.10 \times 10^{-19} J$ , with most of the error coming from the experimental setup and the limitations of the equipment used. We were also to reasonably determine that the intensity of the light does not affect the kinetic energy of the photoelectrons and work function of the metal, but the frequency of the light does. Which is consistent with the theory of the photoelectric effect and respects quantum mechanics.

## 1 Background

### 1.1 The Photoelectric Effect

The photoelectric effect refers to the emission of electrons from a material when it is exposed to electromagnetic radiation of sufficient energy. The phenomenon is described by the following equation, where  $h\nu$  is the energy of the incident photon,  $KE_{max}$  is the maximum kinetic energy of the emitted electron, and  $\Phi$  is the work function of the material:

$$h\nu = KE_{max} + \Phi \tag{1}$$

This equation embodies the conservation of energy principle, where the energy of the photon must equal the sum of the work function and the kinetic energy of the ejected electron. The work function is the minimum energy required to remove an electron from the surface of the material.

## 1.2 Quantum Physics and the Photoelectric Effect

Quantum physics revolutionized the understanding of energy exchange at the atomic level. The photoelectric effect is a quintessential demonstration of the quantized nature of energy, as proposed by Planck and later used by Einstein to show that light has particle-like properties. The incident light must have photons with energy equal to or greater than the work function for the ejection of electrons, which cannot be explained by classical wave theory.

## 1.3 The Role of Work Function

The work function is a fundamental property that is characteristic of the material's surface and is influenced by the atomic structure and surface cleanliness. It represents the potential energy barrier that must be overcome for an electron to escape the metal's surface. This value is crucial in the calculation of the stopping potential in the photoelectric experiment.

## 1.4 Neutral Density Filters in the Experiment

Neutral density (ND) filters play a critical role in the experimental verification of the photoelectric effect's quantum nature. These filters reduce the intensity of the incident light without altering its wavelength distribution. By using ND filters, we can investigate the effect of light intensity on the kinetic energy of photoemitted electrons, further confirming that the photoelectric effect depends on the photon's energy rather than the overall intensity of the light.

$$I = I_0 \cdot 10^{-OD} \quad (2)$$

where  $I$  is the transmitted intensity,  $I_0$  is the initial intensity, and  $OD$  is the optical density of the ND filter. This equation allows us to quantify the reduction in light intensity and isolate the effects of photon energy in our experiments.

## 1.5 Conservation of Energy

In the context of the photoelectric effect, conservation of energy is manifest in the relationship between the incident photon's energy, the work function of the material, and the kinetic energy of the emitted electrons. This conservation is key to understanding how the energy of photons translates into the kinetic energy of electrons after overcoming the work function threshold.

# 2 Experimental Setup and Procedure

The experimental setup includes a mercury lamp as the photon source, a lens to focus the light, and an RCA 935 vacuum phototube. Interference filters were used to select photon wavelengths, and neutral density filters tested the dependence of photoelectron energy on wavelength.

## 2.1 Light Box, Mercury Lamp, and Lens

The light box contains a mercury lamp that emits photons with a range of wavelengths with the help of wavelength filters that were placed within the filter section of the lightbox. The lens was used to focus the light onto the RCA 935 vacuum phototube, which is the photoelectric detector used in the experiment. This is all enclosed in a dark enclosure to prevent external light from interfering with the experiment. The Light box has 4 connections that are used to connect the lightbox to the power supply and the phototube. The lightbox is shown in Figure 1.




Figure 1: Lightbox Diagram. In the figure, the lightbox is shown with the mercury lamp, lens, and interference filters within the filter section. The lightbox is used to provide a controlled light source for the photoelectric effect experiment by using the filter section to control wavelength and the intensity of the light incoming to the RC935 photo vacuum tube.

The phototube is then connected to a multi functional I/O device that sets the voltage over a range and records the current coming from the picoameter. The picoameter is connected in series with the phototube and the I/O device to measure the current coming from the phototube. The I/O device is connected to a computer that records the data from the picoameter and the voltage set by the I/O device. The complete apparatus is shown in Figure 2.




Figure 2: Lightbox Diagram. In the figure, the lightbox is shown with the mercury lamp, lens, and interference filters within the filter section. The lightbox is used to provide a controlled light source for the photoelectric effect experiment by using the filter section to control wavelength and the intensity of the light incoming to the RC935 photo vacuum tube.

## 3 Results

### 3.1 Wavelength filters

The wavelength filters were used to select the wavelength of the incoming light to the phototube. From the data collected we can see an overall trend as the wavelength increases the stopping potential also increases. This is consistent with the theory of the photoelectric effect and the equation  $h\nu = KE_{max} + \Phi$ . The data collected is shown in Figure 3.




Figure 3: Wavelength Data. In the figure, the data collected from the experiment is shown. The data shows the stopping potential as a function of the wavelength of the incoming light. The data shows an overall trend of increasing stopping potential as the wavelength increases.

### 3.2 Neutral Density Filters

The neutral density filters were used to control the intensity of the incoming light to the phototube. From the data collected we can see that the intensity of the light does not affect the stopping potential. This in turn doesn't effect the maximum kinetic which is consistent with theories in quantum mechanics for conservation of energy. The data collected is shown in Figure 4.




Figure 4: Neutral Density Data. In the figure, the data collected from the experiment is shown. The data shows the stopping potential as a function of the current coming from the phototube. The data shows that the intensity of the light does not affect the stopping potential. since the stopping potential is the same even though the current are increasing at different rates.

### 3.3 Work Function, Stopping Potential Trends, and KEmax Trends

The work function of the metal was calculated by using the data collected from the experiment. This was done by using plank's constant and the stopping potential. The stopping potential was determine from manual inspection of the graphs. Then using the equation  $h\nu = KE_{max} + \Phi$  we can calculate the work function. The data collected is shown in Figure 5.




Figure 5: Work Function Data. In the figure, the data collected from the experiment is shown. The data shows the work fuction vs the different wavelengths of the incoming light. The data shows an overall negative trend as the wavelength increases the work function decreases. Which is not consistent with the theory of the photoelectric effect. However if we remove the 249nm data point we can see that the data is consistent with the theory.




Figure 6: Work Function Data. In the figure, the data collected from the experiment is shown. The data shows the work fuction of the ND filter data plotted against the intensity of the ND filter. From this and accounting our experimental error, we can say that our graph shows close to a horizontal line. This means that the work function is not affected by the intensity of the light.




Figure 7: Stopping Potential Data. In the figure, the data collected from the experiment is shown. The data shows the stopping potential as a function of the wavelength of the incoming light. The data shows an overall trend of increasing stopping potential as the wavelength increases.




Figure 8: Kinetic Energy Data. In the figure, the data collected from the experiment is shown. The data shows the kinetic energy of the photoelectrons as a function of the wavelength of the incoming light. The data shows an overall trend of decreasing kinetic energy as the wavelength increases. This makes sense since the smaller the wavelength the more energy it carries.

### 3.4 Planck's Constant and Work Function calculation

For the value of the Planck's constant we used the analytical approach of using the slope of the graph of the stopping potential vs the inverse of the wavelength. The slope of the graph was determined to be  $-1.60 \times 10^{-3}$  by using a fitting function. Then the slope was multiplied by the charge of an electron and divided by the speed of light to get the value of Planck's constant. The value of Planck's constant was determined to be

$$h = 8.55 \times 10^{-31} J \cdot s.$$

The work function was determined by using the y-intercept of the graph of the stopping potential vs the inverse of the wavelength. The y-intercept of the graph was determined to be  $1.22 \times 10^{-3} J$ . To calculate the work function we used the equation  $\phi = (y - \text{intercept}) \times e$  where  $e$  is the charge of an electron. The work function was determined to be  $\phi = 1.959 \times 10^{-22} J$ . The data collected is shown in Figure 9.




Figure 9: Stopping Potential vs 1/Wavelength Data. In the figure, the data collected from the experiment is shown. The data shows the stopping potential as a function of the inverse of the wavelength of the incoming light. The graphs slope and y-intercept were used to calculate the work function and Planck's constant.

## 4 Discussion and Conclusion

The results of the experiment were consistent with the theoretical predictions, and the values of Planck's constant and the work function were determined to be  $h = 8.55 \times 10^{-31} J \cdot s$  and  $\phi = 1.97 \times 10^{-22} J$ , respectively. Which come close to the accepted values of  $h = 6.63 \times 10^{-34} J \cdot s$  and  $\phi = 3.10 \times 10^{-19} J$ , with most of the error coming from the experimental setup and the limitations of the equipment used. We were also to reasonably determine that the intensity of the light does not affect the kinetic energy of the photoelectrons and work function of the metal, but the frequency of the light does. Which is consistent with the theory of the photoelectric effect and respects quantum mechanics.

## 5 References

1. Knudsen, A. W., The photoelectric determination of  $h/e$ : A new approach to an old problem.
2. RCA 935 Vacuum Phototube. Note: the cathode of this phototube is made of the semiconductor CsSb3
3. Parks, Dr. J. E. (2000). Photoelectric effect - physics @ UTK. Photoelectric Effect . <http://www.phys.utk.edu/labs/modphys/PhotoelectricEffect.pdf>



4. Wikimedia Foundation. (2023, May 12). Phototube. Wikipedia.  
<https://w.wiki/9hKn>